

**South Plains College**  
**EMS Programs**  
**Glucometer**

NAME:		DATE:	
COURSE:		SEMSTER:	INITIAL: <span style="border: 1px solid black; display: inline-block; width: 20px; height: 15px;"></span> RETEST: <span style="border: 1px solid black; display: inline-block; width: 20px; height: 15px;"></span>

<b>Task:</b> Assess blood glucose level with a glucometer.
<b>Performance Observations:</b> The candidate shall be able to correctly administer blood glucose level using a glucometer.
<b>Candidate Directive:</b> "Properly assess blood glucose levels with a glucometer."

No.	Task Steps	Pass	Fail
1.	Explain the procedure to the patient.		
2.	Follow manufactures instructors for operation of glucometer.		
3.	Using aseptic technique cleanse site for blood sample.		
4.	Prick finger and gently squeeze to obtain blood specimen.		
4.	Apply blood to test strip.		
5.	Allow monitor to assess blood glucose level.		
6.	Dispose of sharps and contaminated strip in appropriate container.		
7.	Report/document reading of blood glucose level.		

Absolutely satisfied _____ Yes _____ No _____	_____ Pass _____ Fail
Examiner: <span style="border: 1px solid black; display: inline-block; width: 150px; height: 15px;"></span>	Certification # <span style="border: 1px solid black; display: inline-block; width: 100px; height: 15px;"></span>

Documenting Comments: